

Development of two-beam femtosecond/picosecond CARS for high-fidelity thermometry in flames

Bohlin, Alexis

Publication date

2018

Document Version

Final published version

Citation (APA)

Bohlin, A. (2018). *Development of two-beam femtosecond/picosecond CARS for high-fidelity thermometry in flames*. Abstract from 32nd International Congress on High-Speed Imaging and Photonics, Enschede, Netherlands.

Important note

To cite this publication, please use the final published version (if applicable). Please check the document version above.

Copyright

Other than for strictly personal use, it is not permitted to download, forward or distribute the text or part of it, without the consent of the author(s) and/or copyright holder(s), unless the work is under an open content license such as Creative Commons.


Takedown policy

Please contact us and provide details if you believe this document breaches copyrights. We will remove access to the work immediately and investigate your claim.

Development of two-beam fs/ps CARS for high-fidelity thermometry in flames

Alexis Bohlin, Ph.D.

Faculty of Aerospace Engineering, Delft University of Technology



Acknowledgement:





Funding provided by NWO (AES) - Netherlands Organisation for Scientific Research (Vidi grant)



Advanced optical diagnostics are important tools for quantitative combustion analysis

1. Validation/development of a model require multi-parameter diagnostics

3. Measurement challenges
Spatial / temporal / spectral resolution – data acquisition



2. Parameter determination in reacting flows
(e.g. temperature, flow-field and species)

Why use CARS for flame diagnostics?


- Most accurate technique for **thermometry** in reacting flows (wide range of operational conditions).





Inaccuracy $\sim 2-3\%$
Single shot precision $\sim 4-5\%$

- Nanosecond CARS characteristics:

- Non-intrusive, in-situ probe
- High temporal resolution (~ 10 ns)
- High spatial resolution ($\sim 100 \mu\text{m} \times 100 \mu\text{m} \times 1-2$ mm)



- Vibrational CARS, Rotational CARS



Why use CARS for flame diagnostics?


- Most accurate technique for **thermometry** in reacting flows (wide range of operational conditions).



Inaccuracy $\sim 2-3\%$
Single shot precision $\sim 4-5\%$

- Nanosecond CARS characteristics:

- Non-intrusive, in-situ probe
- High temporal resolution (~ 10 ns)
- High spatial resolution ($\sim 100 \mu\text{m} \times 100 \mu\text{m} \times 1-2$ mm)




- Vibrational CARS, Rotational CARS


- Two-beam femtosecond/picosecond CARS

- Picosecond temporal resolution (Near collision independent - Raman linewidths)
- Improved spatial resolution ($40 \mu\text{m} \times 40 \mu\text{m} \times 0.5$ mm)
- 1D and 2D imaging capabilities

Inaccuracy $< 2-3\%$
Single shot precision $\sim 1\%$




Two-beam femtosecond/picosecond CARS




Simplified generic phase-matching-scheme for CARS signal generation


Energy conservation




Phase-matching (momentum conservation)




Spectroscopy in the time-domain




Direct coherent Raman temperature imaging and wideband chemical detection



- Canonical sooting hydrocarbon flat-flame used to benchmark the new techniques.




Burner design (Michelsen group, Sandia)



CARS imaging of flame-wall interaction


- Temperature contour mapping at side-wall quenching burner



Burner design (Dreizler group, TU Darmstadt)


- Motivation

Flame-wall interaction plays a key role in the formation of pollutants in a combustion chamber, such as UHC and CO.




CARS imaging of flame-wall interaction

- Temperature contour mapping at side-wall quenching burner




Burner design (Dreizler group, TU Darmstadt)





- Motivation

Flame-wall interaction plays a key role in the formation of pollutants in a combustion chamber, such as UHC and CO.



Near-wall ultrabroadband CARS imaging: Measurement of thermochemical states





Simultaneous detection of N_2 , O_2 , H_2 , (CO), CO_2 , and CH_4 is achieved.

The excellent imaging resolution allows for thermochemical states of the thermal boundary layer to be probed to within ~ 40 μm of the interface.

In-situ measurement of pressure broadening coefficients

FWI at enhanced turbulence intensities (Work-in-progress)


Turbulence generating grid (blockage 45%, turbulence level $u' / \bar{u} = 6-7\%$),
V-flame operating in the wrinkled flamelet regime





- Single-shot spatially dependent statistics of the 1D flame-front gradient / thickness / position become possible (improving heat transfer models)

Single-shot hyperspectral CARS in the gas-phase


Temperature imaging




Wideband chemical imaging



Dispersive Fourier Transform for MHz detection of CARS/CSRS signals






Synchronized ps/fs laser system for time-resolved non-linear optical spectroscopy/microscopy

- Femtosecond laser (ultrafast amplifier)
7 mJ/pulse @ ~780-810 nm (~35 fs)
- Picosecond laser (SHBC)
2.0 mJ/pulse @ 400 nm (~10 ps)

Distributed auto-ignition combustion modes with reduced NO_x emission



Courtesy of: Arvind Gangoli Rao

Conclusions

- Two-beam femtosecond/picosecond CARS
 - Relevant for 0D, 1D, and 2D temperature measurements in flames when high-fidelity information is needed (inaccuracy <2-3%, precision ~1%)
 - Single-shot quantitative measurements for major species in combustion are within reach (species specific dephasing times, spectroscopy models)
- This ultrafast 1D-CARS technique has been successfully employed at:
 1. Flame-wall interaction burner (head-on and side-wall quenching)
 2. Sooty flames provided on a McKenna burner
- Can this advanced laser diagnostics technique be employed for measurements in engines?
 - Technical challenges for the stability of operation (facility temperature and humidity control, propagating TL-beams through optical ports)